


A Firm Grip with
Precise Operation




Medical Grade Tweezers

Shark Skin Tweezers
“SAMEHADA SESSI”




鮫肌
鑷子

The shark skin processing on the tip and handle enables precise, accurate gripping operations on soft tissues.



Tip



Handle

Merits of Shark Skin Processing

- Enables precise, accurate operations on tissues of lungs, intestinal tracts, membranes, and vascular sheaths.
- Prevents degrading of operational performance resulting from bloodstains on the tip.
- Retains operational performance when surgical gloves are on.
- The weight is engineered with long surgeries in mind and is ideal in such situations.

Order Information

Product Number	Product Name	Specifications	JAN Code
S32222405B0002	Shark Skin Tweezers <line>	Length 240mm Tip Width 2mm	4580469475923
S32222405B0004	Shark Skin Tweezers <thin tipped>	Length 240mm Tip Width 1.5mm	4580469480491
S32223005B0001	Shark Skin Tweezers <line>30cm	Length 300mm Tip Width 3mm	4580469482334

This product was created through the "Ikou No Kizuna" ("Medical-engineering collaboration"), a training course for hybrid medical professionals. (Patent No. 6310051)

Generic Name: Tweezers Medical Equipment Manufacturing and Selling Registration Number: 13B1X0027400038
 ※This is not comply with European Legislation

Manufacturing and Sales Business



株式会社 田中医科器械製作所
 Tanaka Medical Instruments Co., Ltd.

Distributor



株式会社 アステム

Higashiie Business Center, 3-1-2, Higashiie,
 Hakata-ku, Fukuoka City, Fukuoka, 812-0007, Japan
 TEL: +81-924-09-1624

